Caution: Federal law restricts this device to sale by or on the order of a dentist.

THE CARIES DETECTIVE

Caries Detecting Dye

Technical Data

Dark blue dye in propylene glycol liquid formulation

Shelf life: 2 years

Store at room temperature

Avoid direct heat or sunlight

Product Description

Snoop is a propylene glycol based caries detector used to assist the dentist in identifying and removing infected carious dentin. Its dark blue color contrasts with dentin and is easily distinguished from the pulp.

Researchers have reported that carious dentin is made up of two distinct layers. There is an outer layer of infected carious dentin which is soft, discolored, nonvital, nonsensitive and cannot remineralize. This outer layer is stained by Snoop and should be removed. There is also an inner layer of uninfected carious dentin (also referred to as affected dentin) which is vital and sensitive, not infused with bacteria, and may be somewhat softened and demineralized, but which is capable of remineralizing. This inner layer is not stained by Snoop and should not be removed.

Research indicates that the dye in the caries detector bonds to the denatured collagen which is present in the outer infected dentin but which is not present in inner uninfected dentin and normal dentin. Only the outer infected dentin is stained. As a result, Snoop makes it easy to distinguish between outer and inner dentin and enables the dentist to perform a conservative cavity preparation removing only the outer infected dentin.

Warning

- Snoop contains dye that stains. Handle with care.
- Apply Snoop only to the tooth to be stained and not to adjacent teeth.
- If Snoop comes in contact with clothes, carpets or furniture, wash immediately with soap and water.
- If Snoop comes in contact with eyes, flush with copious amounts of water for 15 minutes and seek medical attention.
- Wear protective eye wear. Use rubber dam. Drape patient with plastic bib.
- Keep out of reach of children.

Directions for Use

1. Use of rubber dam is recommended.
2. Use Snoop as early in the cavity preparation procedure as possible to avoid removing too much tooth structure.
3. Rinse treatment area with water and dry.
4. Place applicator straw firmly into dropper spout. Using gentle pressure, carefully apply SNOOP directly to the cavity for 10 seconds.
5. Alternatively, hold bottle vertically and dispense a drop of SNOOP into a mixing well. Using a small sponge, brush or paper point, carefully apply SNOOP to the cavity for 10 seconds.
6. After 10 seconds, rinse the cavity with water while applying suction, then dry.
7. Remove the stained dentin using either a low speed handpiece with a round bur or an excavator.
8. Repeat procedure until no staining occurs.

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PULPDENT Corporation