Instructions For Use

CE
0459

PRODUCT DESCRIPTION

TempCanal Enhanced is a root canal treatment paste recommended for routine use between office visits to cleanse canals and prevent flare-ups, and for use in complicated cases when extended calcium hydroxide treatment is indicated.

The paste flows through 27-gauge x 1" (0.4 x 25mm), 2-side-vent, endodontic irrigation needles for direct, one-handed delivery into the root canal. Closed-end needles prevent dispensing beyond the apex. Non-drying formula will not clog 27-gauge needles.

TempCanal Enhanced has a pH>12 and provides the proven benefits of calcium hydroxide for successful root canal therapy. Use routinely as an intracanal dressing and for treatment of periapical lesions and abscesses, internal and external root resorption, traumatic injuries, root fractures and perforations, apexification, and weeping canals.

INSTRUCTIONS FOR USE

1. Prepare and cleanse canal in the usual manner.
2. Attach luer lock dispenser needle securely onto the syringe. (Alternatively, paste can be placed in the canal using a lentulo paste filler.)
3. Insert needle into the canal to the desired depth. The needle is closed end with two side vents.
4. Depress syringe plunger to dispense paste.
5. Fill canal while slowly withdrawing the needle.
6. Seal access opening with your preferred cement.
7. Discard needle and recap syringe immediately after use
8. TempCanal Enhanced is easily removed with file and irrigation.

Precautions: May be irritating to eyes, skin, respiratory system, gastrointestinal tract, or mucous membranes, especially in large quantities or with prolonged exposure. For use by dental professionals only. Wear gloves and safety glasses. Wash hands after use. Avoid unnecessary exposure.


PULPDENT Corporation

80 Oakland Street • Watertown, MA 02472 • U.S.A.
(617) 926-6666 / (800) 343-4342 / Fax (617) 926-6262
www.pulpdent.com • pulpdent@pulpdent.com • www.activabioactive.com

XP-TE-IN-02
REV: 09/2016

CLINICAL SITUATIONS WHERE CALCIUM HYDROXIDE CAN BE USED

- 1. exudation control
- 2. large periapical lesions
- 3. antibacterial dressing
- 4. temporary root filling
- 5. apical inflammatory resorption
- 6. inflammatory resorption following trauma

- 7. apical internal resorption
- 8. internal-external root resorption
- 9. perforations
- 10. transverse root fractures
- 11. incompletely developed pulpless teeth

(Heithersay, G.S. Calcium Hydroxide in the Treatment of Pulpless Teeth with Associated Pathology, J. Brit Endo Society, 8[2]: 74-93 July 1975)

Note: Apply disposable barrier sleeves/wraps over multiple-use dental dispensers before use with each patient.

For additional information, refer to the following websites:

<http://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/DentalProducts/ucm404472.htm>

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1.0 Commercial Product Name and Supplier

1.1	Commercial product name / designation	TempCanal Enhanced
1.2	Application / Use	Dental material for use by dental professionals.
1.2.2	SIC	851 Human health activity
1.2.3	Use Category	55
1.3	Manufacturer	Pulpdent Corporation 80 Oakland Street, P.O. Box 780 Watertown, MA 02472 USA Telephone: 1 617 926-6666 Fax: 1 617 926-6262 Email: Pulpdent@pulpdent.com
1.4	Emergency Telephone Number	1-800-535-5053 (24 Hour / USA)
1.5	Authorized European Representative	Advena Ltd. Pure Offices, Plato Close Warwick, CV34 6WE United Kingdom

2.0 Hazards Identification

2.1 Classification

2.1.1 Classification according to Regulation (EC) No. 1272/2008 [CLP]

Hazard Class	Hazard Category	Hazard Statement
Skin irritation	2	H315
Eye irritation	2	H319

2.1.2 Classification according to Directive 67/548/EEC

2.1.3 Classification according to WHMIS

2.2 GHS Label Elements

Hazard Pictograms

Signal Word: **WARNING**

Restricted to use by dental professional only.

Hazard Statements H319: Eye irritation. 2. May cause eye irritation.

H315: Skin irritation. 2. May cause skin irritation.

Precautionary Statements P280: Wear protective gloves and eye protection

P305+P351+P338: If IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If in eyes, rinse cautiously with water for several minutes.

P337+P313: If eye irritation persists, get medical advice/attention.

P302+P352: If on skin, wash with plenty of soap and water.

3.0 Composition

3.1 Chemical characterization of the preparation Calcium hydroxide suspended in a cellulose/water matrix.

3.2	Hazardous ingredients				
	CAS Number	Name of the Ingredient	Concentration	Classification per 67/548/EEC	Classification per Regulation (EC) No.1272/2008 (CLP)
	1305-62-0	Calcium Hydroxide	25%	Irritant (Xi) R 36 / 38	Skin irritation; 2 Eye irritation; 2
4.0	First Aid Measures				
4.1	Special Instructions	Show this safety data sheet to medical personnel. Get medical attention in case of uncertainty. This product presents minimal hazard in the quantity of material provided and typically used. Only with prolonged exposure is it irritating to mucous membranes, eyes, skin, gastrointestinal tract or respiratory tract.			
4.2	Eye Contact	Keep eyelids apart; flush with running water for 15+ minutes. Get medical attention for persistent irritation.			
4.3	Skin Contact	Remove any contaminated clothing. Wipe off excess from skin. Immediately flush skin water. Get medical attention if irritation persists.			
4.4	Ingestion	Rinse mouth with water. Do not induce vomiting. Give large amounts of water to dilute. Get immediate medical attention for ingestion of large amount. Never give anything by mouth to an unconscious person.			
4.5	Inhalation	Move to fresh air. Administer artificial respiration as needed; seek medical attention.			
4.6	Precautions for first responders	Wear safety glasses and gloves to prevent contact.			
4.7	Information for physician				
	Symptoms	Irritation, pain or redness in eyes, on skin, on mucous membranes.			
	Hazards	May be irritating to eyes, skin, gastrointestinal tract or mucous membranes especially with prolonged exposure There is a risk of serious damage to eyes.			
	Treatment	Same as First Aid above.			
5.0	Fire Fighting Measures				
5.1	Suitable extinguishing media	Not a fire hazard. Extinguish fire with agent suitable for surrounding fire.			
5.2	Extinguishing media to avoid	None			
5.3	Special exposure hazards in a fire	None			
5.4	Special protective equipment: fire-fighters	Self-contained breathing apparatus			
6.0	Accidental Release Measures				
6.1	Personal precautions.	Wear gloves, lab coat and chemical splash goggles.			
6.2	Environmental precautions	None			
6.3	Method for clean up	Wipe up or absorb with paper towel. Wash area of spill with water-saturated paper towels or cloths. Place spilled material and absorbent towels in container for disposal.			
7.0	Handling and Storage				
7.1	Handling	For use by dental professionals only. Wear gloves, safety glasses; wash hands after use. Avoid unnecessary exposure. Follow good hygiene practices. Do not smoke, eat or drink. Wash hands after use.			
7.2	Storage	Keep tightly capped in original container. Store at cool room temperature. Avoid extremes of temperature			

		(>27°C/80°F, <5°C/40°F) and acids. Dental preparation used for cavity lining
7.3	Specific uses	

8.0 Exposure Controls / Personal Protection

8.1	Exposure limit values	Not applicable
8.2	Exposure controls	
8.2.1	Occupational exposure controls	Wear safety glasses, lab coat and gloves. No other equipment required under normal conditions of use of this product in the quantity provided.
8.2.1.1	Respiratory protection	Good general ventilation is sufficient to control any airborne vapors.
8.2.1.2	Hand protection	No special requirements other than the usual surgical gloves.
8.2.1.3	Eye protection	No special requirements other than the usual safety glasses.
8.2.1.4	Skin protection	No special requirements. Good personal hygiene and safety practices and wearing a lab coat should protect dental staff from unnecessary exposure.
8.2.1.5	Other controls	Emergency eye wash fountain should be close by. Wash hands after use.
8.2.2	Environmental exposure controls	Avoid release into environment; calcium hydroxide may cause pH variation.

9.0 Physical and Chemical Properties

9.1	Characteristics	
9.1.1	Appearance /Color / Physical state	Thin white paste
9.1.2	Odor	Mild, characteristic
9.2	Important health, safety and environmental information	
9.2.1	pH	pH > 12
9.2.2	Boiling point	Decomposes
9.2.3	Flash point	Not combustible
9.2.4	Flammability (solid, gas)	Not combustible
9.2.5	Explosive properties	Not applicable
9.2.6	Oxidizing properties	Not determined
9.2.7	Vapor pressure	< 1 mm Hg (Id: A)
9.2.8	Specific gravity	1.470 g/ml
9.2.9	Solubility in water	Slight
9.2.10	Partition coefficient	Not determined
9.2.11	Viscosity	Not determined
9.2.12	Vapor density	Not determined
9.2.13	Evaporation rate	Not determined

10.0 Stability and reactivity

10.1	Conditions to avoid	Extremes of temperature (>27°C/80°F, <5°C/40°F).
10.2	Materials to avoid	Calcium hydroxide reacts violently with maleic anhydride, nitro-ethane, nitro-methane, nitro-paraffin, nitro-propane, phosphorus. As a strong alkali, it is incompatible with acids.
10.3	Hazardous decomposition products	Caustic fumes of calcium oxide when heated to decomposition.

10.4	Further information	Stable under normal conditions.
11.0 Toxicological information		
11.1	Acute toxicity	Not toxic. Minimal health hazard under normal conditions of use.
11.2	Irritation and corrosiveness	Minimal health hazard under normal conditions of dental practice. Prolonged exposure may cause irritation of mucous membranes, eyes, skin, gastrointestinal tract or respiratory tract.
11.3	Sensitization	Not applicable.
11.4	Sub-acute, sub-chronic and prolonged toxicity	Narrowing of the esophagus may occur weeks to years after ingestion of large quantities, making swallowing difficult. Persons with pre-existing skin conditions may be more susceptible to the irritating effects of calcium hydroxide
11.5	Carcinogenicity, Mutagenicity, Reproductive Toxicity	Not considered a carcinogen, mutagen or teratogen. Not toxic to the reproductive system.
11.6	Empirical data	Not available
11.7	Clinical Experience	Pulpdent calcium hydroxide preparations have been used safely and effectively by dental professionals since 1947.
12.0 Ecological Information		
12.1	Ecotoxicity	Not considered detrimental to water or soil. Avoid release into environment as calcium hydroxide may cause pH variation. Follow government regulations.
13.0 Disposal Considerations		
13.1	Regulations	Follow all local and national government regulations in disposing material or contaminated packaging.
14.0 Transport Information		
14.1	Restrictions	None. Not regulated.
14.2	UN Number	None
14.3	Technical name	Not applicable
14.4	Packing group	Not applicable
14.5	IATA class	Not applicable
15.0 Regulatory Information		
15.1	EU Regulations	Class IIa medical device under MDD 93/42/EEC.
15.2	US FDA	Class II medical device
15.3	Health Canada	Class III medical device
16.0 Other information		
16.1	List of relevant R phrases	R36/38: Irritating to eyes and skin.
16.2	Hazard Statements	H319: Eye irritation. 2. May cause eye irritation. H315: Skin irritation. 2. May cause skin irritation.
16.3	Precautionary Statements	P280: Wear protective gloves and eye protection P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If in eyes, rinse cautiously with water for several minutes. P337+P313: If eye irritation persists, get medical advice/attention.

	P302+P352: IF ON SKIN, wash with plenty of soap and water. TempCanal Enhanced is to be sold to and used by dental professionals only.
16.4	Restrictions on use
16.5	Further information The information presented herein is believed to be factual as it has been derived from the works of persons believed to be qualified experts. However, nothing contained in this information is to be taken as a warranty or representation for which Pulpdent Corporation bears legal responsibility. The user should review any recommendations in the specific context of the intended use to determine whether they are appropriate.
16.6	Sources of key data National Institute for Occupational Safety (NIOSH) Occupational Safety and Health Administration (OSHA) Eur-Lex European Union Law: Regulation (EC) No. 1272/2008 (CLP) and Regulation (EC) No. 1907/2006 (REACH). Guidance on the compilation of safety data sheets. Version 1.1; December 2011. European Chemicals Agency
16.7	Information which has been added, deleted or revised. This Safety Data Sheet has been revised to meet the requirements of the GHS SDS format, Regulations (EC) No. 1272/2008 (CLP) and (EC) No. 1907/2006 (REACH). Specifically, Sections 2.1, 2.2, 3.2, 16.2, 16.3 have been modified.

PULPDENT Corporation

80 Oakland Street • Watertown, MA 02472 • U.S.A.
(617) 926-6666 / (800) 343-4342 / Fax (617) 926-6262
www.pulpdent.com • pulpdent@pulpdent.com • www.activabioactive.com