# Lime-Lite™

## Light Cure Cavity Liner

- Fluoride Releasing
- Radiopaque
- Dentin Shade
- Contains Hydroxyapatite

**CONTENTS:**

- 4 1.2ml syringes Lime-Lite™  
8 Dispenser tips

Instructions for use


MSDS

Caution: Federal law restricts this device to sale by or on the order of a dentist.

Store in a cool, dry place

MADE IN USA

For use with all dental restorative materials and cements

**Lime-Lite™**  
Light Cure Cavity Liner

**PULPDENT®**

# Lime-Lite™ Light Cure Cavity Liner

- Fluoride Releasing • Radiopaque • Universal Shade • Contains Hydroxyapatite

Lime-Lite™ is a light cure, fluoride releasing, radiopaque cavity liner and base material specially formulated for use with adhesives, composites and conventional restorative materials. Lime-Lite™ contains hydroxyapatite in a urethane dimethacrylate resin. It chemically bonds to adhesive primers, composites and other resin based materials, and it micromechanically bonds to dentin.

Lime-Lite™ releases favorable Calcium ions, Hydroxyl ions, Fluoride ions and Phosphate ions. These are known to be beneficial to tooth structure, to stimulate secondary dentin formation and to have cariostatic properties. Lime-Lite™ sets very hard upon light curing, and it is virtually insoluble in water and oral fluids.

Lime-Lite™ can be applied directly from the syringe through pre-bent dispenser tips for precise placement or can be dispensed on a pad or mixing well and placed on the tooth with a suitable instrument.

**TECHNICAL DATA**

Active ingredients: calcium hydroxy phosphate (hydroxyapatite), urethane dimethacrylate resin, barium sulfate, fluoride salt, photo initiator.

- Visible Light Cure
- Universal Shade
- Store at room temperature
- Do not freeze
- Shelf Life: 2 years

**WARNING:** Wash uncured Lime-Lite™ from skin and soft tissue. If Lime-Lite™ contacts eyes, wash with copious amounts of water. Seek medical attention, if necessary.

**INSTRUCTIONS FOR USE**

1. Prepare cavity, rinse and dry.
2. Place a pink 20 gauge dispenser tip securely onto Lime-Lite™ syringe.
3. Apply Lime-Lite™ to the dentin in increments. Pull back on syringe plunger to stop flow of material. The first layer should be not more than 1/2 mm thick. Remove Lime-Lite™ from enamel margins and light cure for 20-30 seconds, depending on the light you are using.
4. Apply additional layers, if indicated, in increments of 1 1/2 mm – 2 mm. Remove Lime-Lite™ from enamel margins and light cure each layer for 30 seconds. Cap syringe immediately after use.
5. Proceed with restoration. Lime-Lite™ is compatible with all dental restorative materials. **NOTE:** Lime-Lite™ has a urethane dimethacrylate base and is similar in chemistry to dental composites. It chemically bonds to composites and adhesive primers.

When increased bond strength to dentin is desirable, apply DenTASTIC or your preferred dental adhesive to the dentin in accordance with the instructions for use for the adhesive primer. Then apply Lime-Lite™ over the adhesive primer following instructions # 2-5 above.

For direct pulp capping, we recommend placing Pulpdent Paste or Multi-Cal on the exposed pulp. Lime-Lite™ can be applied over the pulpal dressing and light cured.

If using Lime-Lite™ under water based cements, such as zinc phosphate, polyacrylic or glass ionomer cement, remove any uncured resin from the Lime-Lite™ surface (oxygen inhibited layer) by wiping it with a pected of cotton moistened with alcohol.

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## PULPDENT Corporation