# PULPDENT

## Spee-Dee™ Build-Up

### Multi-Purpose Core and Tooth Build-Up Resin

Dual Cure • Fluoride Releasing • Hydrophilic  
Radiopaque • Self-adhesive • Self-etching to Dentin

Spee-Dee Build-Up is a multi-purpose, dual cure, resin-based dental restorative and core build-up material. Formulated with Pulpdent's proprietary Embrace resin chemistry, it is moisture tolerant and has self-adhesive and self-etching properties that are compatible with all bonding agents. The material is fluoride releasing and radiopaque. It contains no Bisphenol A, no Bis-GMA, and no BPA derivatives.

Spee-Dee Build-Up is dense like dentin and therefore cuts like dentin. Spee-Dee's acidic monomer interacts with the urethane resin to produce a greater conversion of double bonds and a more complete cure, resulting in a denser material with increased hardness. As a result, rotary instruments transition smoothly from tooth to Spee-Dee, cut accurately, and do not ditch, gouge, gum up or chatter.

Spee-Dee Build-Up is used for core build-up, post cementation, and as a multi-purpose repair material for fractured cusps and broken down teeth, defects in crown preparations, lost partial or whole fillings, and chipped or broken crowns or bridges.

Posts and cores can be placed in one step ensuring a homogeneous, one-piece internal structure. Spee-Dee Build-Up also serves as an ideal core material in those cases where no post is indicated (as in fractures) because it flows into cracks and voids, bonds to tooth structure and is fluoride-releasing.

In restorations with retention form, Spee-Dee Build-Up's self-adhesive and self-etching properties may be sufficient for post and core procedures without the use of a bonding agent. However, to achieve highest bond strengths, etching and bonding agents are recommended.

#### Physical Properties

Compressive strength: 44,080 psi / 304 MPa

Flexural strength: 14,065 psi / 97 MPa

Self-cure intra-oral setting time at 37°C: 3:15 minutes

Light cure setting time: 20 seconds for a halogen light with a minimum of 300 mW/cm<sup>2</sup>. Depth of cure = 4 mm. Cures with all lights. More powerful lights will cure faster.

#### How to use the double barrel cartridge:

1. Remove cap. If necessary bleed the cartridge so that base and catalyst are at the orifice of the barrels. Place a mixing tip on the double barrel cartridge.
2. To ensure an even mix of base and catalyst, dispense 2-3 mm of material onto a pad and discard this material.
3. Dispense evenly mixed material directly onto the tooth or into the restoration.
4. Discard mixing tip. Recap cartridge. Do not cross-contaminate base and catalyst.

#### DIRECTIONS FOR USE:

To achieve highest bond strengths, etching and bonding agents are recommended. If Spee-Dee Build-Up is used without etching and bonding agents, lightly dry tooth surface and remove excess water with compressed air or a cotton pellet, and apply Spee-Dee Build-Up to the slightly moist tooth. Slightly moist tooth surfaces exhibit neither dryness nor pooling of water. Tooth surfaces should be shiny or glossy. Overly wet tooth surfaces will result in decreased bond strength. Ceramic, metal, resin desensitized, prehybridized and cured resin surfaces should be dry.

#### Vital Tooth Build-up

1. Remove the old restoration and/or decay.
2. To achieve highest bond strengths, etch and place bonding agent according to the manufacturer's instructions. Light cure. A dual cure bonding agent is recommended with self-cure and dual cure resin core materials.
3. Remove cap from double barrel cartridge. If necessary bleed the cartridge so that base and catalyst are at the orifice of the barrels.
4. Place a mixing tip on the double barrel cartridge and insert into an appropriate 1:1 cartridge dispenser. To ensure an even mix of base and catalyst, dispense 2-3 mm of material onto a pad and discard this material.
5. Dispense material into the preparation site. Light cure for 20 seconds. The material will also self-cure. Carefully recap cartridge. Do not cross-contaminate base and catalyst.
6. Prepare build-up to receive the restoration.

#### Non-Vital Tooth Post and Core Build-up

1. Prepare post space using post system of choice.
2. To achieve highest bond strengths, etch and place dual cure bonding agent according to manufacturer's instructions.
3. Remove cap from double barrel cartridge. If necessary bleed the cartridge so that base and catalyst are at the orifice of the barrels.
4. Place a mixing tip on the double barrel cartridge and insert into an appropriate 1:1 cartridge dispenser. To ensure an even mix of base and catalyst, dispense 2-3 mm of material onto a pad and discard this material.
5. Dispense Spee-Dee Build-Up into the post hole and place the post. Light cure for 20 seconds. The material will also self-cure.

6. Build-up the core using “Crown Up” technique by placing material, curing, adding more material, curing, and repeating this until the desired build-up is accomplished. Carefully recap cartridge. Do not cross-contaminate base and catalyst.
7. Prepare the build-up to receive the restoration.

## For Repairs

### Temporary Patch Material

1. Clean the area to be patched (anesthesia usually not necessary).
2. Etch and place bonding agent in accordance with manufacturer’s instructions.
3. Place Spee-Dee Build-Up as a temporary patch and light cure 20 seconds.
4. Check and adjust occlusion, if necessary.
5. Reschedule patient for appropriate permanent restoration.

### Recurring Decay Under an Existing Crown Margin

1. Remove decay.
2. Etch and bond internal surface of cavity preparation.
3. Place Spee-Dee Build-Up as a restorative material and light cure, finish and polish.

### Broken Temporary Bridge

1. Clean area to be repaired.
2. Using a carbide bur, cut a preparation in the fracture area of the temporary bridge creating a space for Spee-Dee Build-Up.
3. Place Spee-Dee Build-Up and light cure.

### Warning:

- For professional use by dental personnel only.
- Pulpdent Spee-Dee Build-Up may cause irritation to eyes and skin on contact. Take all necessary precautions to avoid contact with eyes and soft tissue. If materials contact eyes, rinse with copious amount of water. If materials contact skin, wash skin with soap and water. Seek medical attention if irritation persists. See MSDS for further information.
- Wear protective eyewear and gloves. Wash hands after use. Emergency eye wash fountain should be available.

### Note:

Apply disposable barrier sleeves/wraps over multiple-use dental dispensers before use with each patient. For additional information, refer to:

<http://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/DentalProducts/ucm404472.htm>

### Storage

Store tightly capped in original container at cool room temperature.  
Avoid direct sunlight, extremes of temperature and cross-contamination between parts.

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### **PULPDENT Corporation**
80 Oakland Street • Watertown, MA 02472 • U.S.A.